



## SECTION 1: IDENTIFICATION

Product Name: Boule Con-5Diff G2  
Part Number: 502-115

**Application of the substance / mixture:** For In Vitro Diagnostic Use

### Manufacturer/Supplier

Address: Produced for Boule Medical AB by:  
Clinical Diagnostic Solutions, Inc.  
1800 NW 65th Avenue  
Plantation FL, 33313, USA  
www.cdsolinc.com

Phone No: +1 954 791 1773  
Fax No: +1 954 791 7118

## SECTION 2: HAZARD(S) IDENTIFICATION

### Classification of the substance or mixture

#### GHS Classification in accordance with 29 CFR 1910

This mixture does not have any hazardous classification.

The mixture has been reviewed in accordance with 29 CFR 1910. The product does not contain any substances at concentrations considered to be hazardous to human health or the environment.

## SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

CAS-No.	Component	Classification	Concentration
-	Proprietary Mixture	-	-

## SECTION 4: FIRST AID MEASURES

### Description of first aid measures

#### General advice

Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

#### In case of skin contact

Wash off with soap and plenty of water. Consult a physician.

#### In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

#### If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.



**Most important symptoms and effects, both acute and delayed**

None

**Indication of any immediate medical attention and special treatment needed**

No data available

## SECTION 5: FIREFIGHTING MEASURES

· **Suitable extinguishing agents**

Dry chemical, carbon dioxide (CO<sub>2</sub>), water spray or regular foam.

- Caution: CO<sub>2</sub> will displace air in confined spaces and may cause an oxygen-deficient atmosphere.

· **Special hazards arising from the substance or mixture**

There are no unique chemical or reactivity hazards that would impact firefighting decisions due to the chemicals in this product.

## SECTION 6: ACCIDENTAL RELEASE MEASURES

· **Personal precautions, protective equipment and emergency procedures**

Minimize exposure by using appropriate personal protective equipment as listed in Section 8. Stop leak if possible. Keep unprotected people away.

· **Environmental precautions**

Prevent liquid from entering sewage system, storm drains, surface waters, and soil.

· **Methods and material for containment and cleaning up**

Block small volumes of spilled or spattered product with paper towels or similar materials.

- Contain larger spills by placing absorbents around the outside edges of the spill. Absorb with any material Suitable for water-based liquids - e.g. paper towels, universal sorbents, sand, diatomite, sawdust, etc.

**Clean the affected area. Suitable cleaners are:**

- Detergent or similar cleansing agent.

Dispose of spilled and contaminated material in Biohazard Waste. See Section 13 for information that may impact disposal of materials contaminated with this product.

## SECTION 7: HANDLING AND STORAGE

· **Precautions for safe handling:** Use personal protective equipment required by the Laboratory.

· **Measure to prevent fire:** No special measures required.

· **Storage requirements:** Store only in the original container. Store in a cool (2 - 10 °C), away from heat and direct sunlight.

· **Further information about storage conditions:** Protect from heat and direct sunlight.



## SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

### Components with Occupational Exposure Limits

Substance name	Concentration	CAS No	Occupational exposure limits
None listed			

#### Exposure controls

##### Personal protective equipment

Always maintain good housekeeping and follow general precautionary measures. Do not eat, drink or store food and beverages in areas where chemicals or specimens are used. Wash hands before breaks, after handling reagents and specimens, and at the end of the work shift.

##### · Breathing equipment:

Normal use and storage of product - respiratory protection is not necessary if room is well ventilated.

**Small-volume spills** (e.g. small enough to clean up with a paper towel or small sorbent pad) - respiratory protection is not necessary.

##### Hand protection:

Wear water-resistant gloves if hand contact with the material is anticipated. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices.

##### · Material of gloves and breakthrough time of the glove material:

The glove material must be suitable for use in a clinical chemistry laboratory and have a measured breakthrough time of at least 30 minutes, such as those with a Class 2 protection index per EN374 (or equivalent standard applicable in your region). NOTE: This recommendation applies only to the product stated in this Safety Data Sheet. When dissolving in or mixing with other substances, contact the supplier of approved gloves.

##### · Eye protection:

Wear safety glasses or other protective eyewear.

##### · Body protection:

Normal use: protect personal clothing from spatters and small spills. Wear a laboratory coat (or other protective clothing required by your institution). Larger spills (e.g. that can saturate cloth): wear appropriate water-repellant covering over clothing.



## SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

Property	Value
Physical state	Thick Liquid
Color	Red
Odor	Odorless
pH-value at 20 °C (68 °F)	Not applicable
Melting point	Not determined
Boiling point	Not determined
Flash point	Not applicable
Flammability (solid, gas)	Non flammable
Auto igniting	Product is not self-igniting.
Danger of explosion	Product does not present an explosion hazard.
Vapor density	Not applicable
Density at 20 °C (68 °F)	Not determined
Solubility in water	Completely soluble

## SECTION 10: STABILITY AND REACTIVITY

- **Thermal decomposition / conditions to be avoided**  
No decomposition if used and stored according to specifications.
- **Possibility of hazardous reactions:** Reacts with strong acids and bases.
- **Conditions to avoid:** Sunlight, warmth, heat and fire.
- **Incompatible materials:** Strong acids and bases.
- **Hazardous decomposition products:** No dangerous decomposition products known.

## SECTION 11: TOXICOLOGICAL INFORMATION

### Acute toxicity

No data available

### Dermal

No data available

### Skin corrosion/irritation

No data available

### Serious eye damage/eye irritation

No data available

### Respiratory or skin sensitization

No data available

### Germ cell mutagenicity

No data available



## SECTION 12: ECOLOGICAL INFORMATION

- **Toxicity:** No further relevant information available.
- **Marine pollutant:** No
- **Additional Ecological Information**
- **General notes:** Do not allow undiluted product or large quantities of it to reach ground water, water course, or sewage system.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable
- **vPvB:** Not applicable

## SECTION 13: DISPOSAL CONSIDERATIONS

### HANDLING OF WASTE

- General information: Before handling waste, see section 8 Exposure controls/Personal protection.  
This product should be disposed in Biohazard waste or infectious medical waste. Disposal by incineration is recommended.  
Biohazard waste must be transported by an approved transporter. For regular transport of Biohazard waste, the user is responsible for providing a transport document. Dispose in accordance with federal, state and local regulations.
- Handling of packaging: Do not clean or recycle. Refer to applicable local regulations and institutional policies.

## SECTION 14: TRANSPORT INFORMATION

- DOT (US):**  
Not dangerous goods
- IMDG:**  
Not dangerous goods
- IATA:**  
Not dangerous goods



## SECTION 15: REGULATORY INFORMATION

• **Section 302/304 (40CFR355.30 / 40CFR355.40):**

The product does not contain listed substances.

• **Section 313 (40CFR372.65):**

The product does not contain listed substances.

• **California Proposition 65 (USA):**

• **Chemicals known to cause cancer:**

None of the ingredients is listed.

• **Chemicals known to cause female reproductive toxicity:**

None of the ingredients is listed.

• **Chemicals known to cause male reproductive toxicity:**

None of the ingredients is listed.

• **Chemicals known to cause developmental reproductive toxicity:**

WARNING! This product contains as chemical known to the State of California to cause birth defects or other reproductive harm.

Neomycin Tridulfate Salt Hydrate CAS 1405-10-3

WARNING! This product contains as chemical known to the State of California to cause developmental harm.

Cycloheximide CAS 66-81-9

### Pennsylvania Right to Know Components

Procaine Hydrochloride	CAS 51-05-8
Edetate Disodium Dihydrate	CAS 6381-92-6
Sodium [2S-(2 $\alpha$ ,5 $\alpha$ ,6 $\beta$ )]-3,3-dimethyl-7-oxo-6-9phenylacetamido0-4-thia-1-azabicyclo[3.2.0]heptane-2-carboxylate	CAS 69-57-8
Citric Acid	CAS 77-92-9
2-Iodoacetamide	CAS 144-48-9
Sodium hydroxide	CAS 1310-73-2
Sodium dihydrogenorthophosphate monohydrate	CAS 10049-21-5
Polyethyglycol	CAS 9003-11-6
Sodium Chloride	CAS 7647-14-5
N-Carbamoylmethyliminodi	CAS 26239-55-4
Mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1)	CAS 55965-84-9
Dihydrostreptomycin Sulfate	CAS 5490-27-7
Neomycin Tridulfate Salt Hydrate	CAS 1405-10-3
Cycloheximide	CAS 66-81-9

### New Jersey Right to Know Components

Procaine Hydrochloride	CAS 51-05-8
Edetate Disodium Dihydrate	CAS 6381-92-6
Sodium [2S-(2 $\alpha$ ,5 $\alpha$ ,6 $\beta$ )]-3,3-dimethyl-7-oxo-6-9phenylacetamido0-4-thia-1-azabicyclo[3.2.0]heptane-2-carboxylate	CAS 69-57-8
Citric Acid	CAS 77-92-9
2-Iodoacetamide	CAS 144-48-9
Sodium hydroxide	CAS 1310-73-2
Sodium dihydrogenorthophosphate monohydrate	CAS 10049-21-5
Polyethyglycol	CAS 9003-11-6
Sodium Chloride	CAS 7647-14-5
N-Carbamoylmethyliminodi	CAS 26239-55-4
Mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1)	CAS 55965-84-9
Dihydrostreptomycin Sulfate	CAS 5490-27-7
Neomycin Tridulfate Salt Hydrate	CAS 1405-10-3
Cycloheximide	CAS 66-81-9



## Massachusetts Right to Know Components

Sodium hydroxide

CAS 1310-73-2

Cycloheximide

CAS 66-81-9

## SECTION 16: OTHER INFORMATION

### Abbreviations and acronyms

CAS No:	Chemical Abstracts Service number.
CMR properties:	Carcinogenic, Mutagenic or toxic for Reproduction
STOT SE	Specific target organ toxicity - single exposure
LC <sub>50</sub> :	Lethal Concentration. In ecotoxicology, the LC <sub>50</sub> is the concentration which kills 50 % of a population of one species, within a specified period of time.
LD <sub>50</sub> :	Lethal Dose. The LD <sub>50</sub> is the dose of a substance which kills 50 % of a population of one species and is expressed as weight (mg, g) or as weight per weight of test animal (mg/kg).
pH:	pH is a measure of the acidity or basicity of an aqueous solution.
pKa:	The symbol for the acid dissociation constant at logarithmic scale.
ppm:	parts per million.
vPvB substance:	Very persistent and very bioaccumulative substance.
WEL:	Workplace Exposure Limits.

### Contact

General information about this product:

Clinical Diagnostic Solutions, Inc.

1800 NW 65<sup>th</sup> Ave.

Plantation, FL 33313